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1. Introduction

The SEDILAB tester was developed to test the hygienic conditions of milk production in dairy farms.

Underlying is the following observation: the more particles and dirt in milk, the poorer hygienic conditions during milking, pouring and cooling.

2. Principle

Raw milk flows through a filter. The particles in the milk are trapped in the filter. The filter is thereafter dried and the results are compared to reference values that each farm or dairy factory should have prepared. (Funke-Dr. N. Gerber's reference tables are suitable for this purpose and can be ordered under art. no. 4920)

3. Installation

Connect the tube of the rubber ball to the unit (on the top of it, to the small pipe). Fix the unit to a stable support (e.g. table or lab bench) by tightening the vice. The outlet of the discharge tubing should be directed to a waste or suitable waste collector. The filter card should be positioned on top of the sieve.

4. Using Sedilab

Tip up the cylinder to position the filter.

Transfer the milk which has been tested in the milk cylinder. Close the top of the cylinder using the palm of the left hand. Apply pressure by pressing the rubber ball with the right hand. The milk will flow through the filter. Lay the filter on a dry surface for drying. Compare the result to the standard values and classify the sample.

5. Cleaning

Wash the unit with detergent and water at the end of each group of tests. The rubber sealing should be frequently controlled and periodically replaced.

6. Note

Raw milk should not be filtered before testing. Should this happen, the test would not come to any reliable results.